

Michigan ADULT TREATMENT SEPSIS

Initial Date: 5/31/2012 Revised Date: 10/25/2017

Section 3-5

Sepsis

It is the purpose of this policy to recognize and treat sepsis early to promote optimal care and survival of patients who may be septic. This protocol applies to patients 14 years and above with a clinical suspicion of systemic infection who have 2 or more of the inclusion criteria. These patients are defined as meeting criteria for suspicion of sepsis and should be evaluated and treated per this protocol.

INCLUSION CRITERIA

- 1. Clinical suspicion of systemic infection, and two or more of the following:
 - A. Hyperthermia temp >38° C (100.4 F)
 - B. Hypothermia temp<36° C (96.8 F)
 - C. Heart rate >90bpm
 - D. Respiratory rate <10 or >20 perminute
 - E. SBP <90 mmHg or evidence of hypoperfusion

Treatment

- Follow General Pre-Hospital Care protocol.
- 2. Place patient in supine position.
- 3. Start large bore IV catheter.
 - 4. Start second large bore IV catheter, if time permits.
- (5. Place on cardiac monitor and treat rhythm according to appropriate protocol.
 - Place on continuous pulse oximetry.
 - Measure blood glucose.
 - 8. If the patient meets inclusion criteria, administer a NS IV/IO fluid bolus up to 1 liter, wide open. Reassess the patient, repeat boluses to a maximum of 2 L NS as long as vital sign abnormalities persist.
 - 9. If hypotension persists, refer to Shock Protocol.
 - (Optional) Measure ETCO2 level. If CO2 < 25, report level to the receiving facility as soon as possible.

Michigan ADULT TREATMENT SEPSIS


Initial Date: 5/31/2012 Revised Date: 10/25/2017

Section 3-5

Follow General Prehospital Care Protocol

Clinical suspicion of systemic infection and two or more of the following:

- Hyperthermia temp > 38 ° C (100.4° F)
- Hypothermia temp < 36° C (96.8° F)
- Heart rate >90bpm
- Respiratory rate <10 or >20 per minute
- SBP <90 mmHg or evidence of hypoperfusion



MCA Name: Bay County MCA MCA Board Approval Date: 02/26/2018 MCA Implementation Date: 05/01/2018 Protocol Source/References: